


Anti-Qprt Polyclonal Antibody For Wb And Elisa - Q15274

Item Number: CM0017587



Introduction

High-quality rabbit polyclonal antibody against human QPRT (nicotinate-nucleotide pyrophosphorylase). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for quinolinic acid catabolism research.

[Learn More](#)

Parameter	Value
Product Name	QPRT Rabbit pAb
Remarks/Alias	QPRTase; HEL-S-90n; QPRT
Species	Human
GeneID (Human)	23475
GeneID	23475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human QPRT (NP_055113.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15274
Protein Weight	31kDa
Shipping	Ice bag