

Anti-Qsox1 Polyclonal Antibody For Wb, Elisa - 000391

Item Number: CM0028567



Introduction

QSOX1 rabbit polyclonal antibody validated for Western blot and ELISA with human and mouse reactivity. Targets sulfhydryl oxidase 1, a key enzyme in disulfide bond formation and extracellular matrix remodeling.

[Learn More](#)

Parameter	Value
Product Name	QSOX1 Rabbit pAb
Remarks/Alias	Q6; QSCN6; QSOX1
Species	Human
Gene ID	5768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human QSOX1 (NP_001004128.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00391
Protein Weight	83kDa
Shipping	Ice bag