

Anti-Qsox2 Polyclonal Antibody For Wb And Elisa - Q6Zrp7

Item Number: CM0018700

Introduction

QSOX2 Rabbit pAb detects human QSOX2, cross-reacts with mouse, and is validated for WB and ELISA. Ideal for studying sulfhydryl oxidase 2 in apoptosis and disulfide bond formation.

[Learn More](#)

Parameter	Value
Product Name	QSOX2 Rabbit pAb
Remarks/Alias	SOXN; QSCN6L1; QSOX2
Species	Human
GeneID (Human)	169714
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 370-650 of human QSOX2 (NP_859052.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q6ZRP7
Protein Weight	78kDa
Shipping	Ice bag