

Anti-R9Ap/Rgs9Bp Polyclonal Antibody For Wb And Elisa - Q6Zs82

Item Number: CM0019204

Introduction

Rabbit polyclonal antibody targeting R9AP/RGS9BP, a membrane anchor for RGS9-1 in phototransduction. Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Ideal for vision and GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	R9AP/RGS9BP Rabbit pAb
Remarks/Alias	R9AP; RGS9; PERRS; R9AP/RGS9BP
Species	Human
Gene ID (Human)	388531
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 1-107 of human R9AP/RGS9BP (NP_997274.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZS82
Protein Weight	25kDa
Shipping	Ice bag