

Anti-R9Ap/Rgs9Bp Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q6Zs82

Item Number: CM0004047

Introduction

Rabbit monoclonal antibody targeting human R9AP/RGS9BP (Q6ZS82). Validated for WB, IF/ICC, ELISA. Critical for studying visual phototransduction and GPCR signaling regulation. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	R9AP/RGS9BP Rabbit mAb
Remarks/Alias	R9AP; RGS9; PERRS; R9AP/RGS9BP
Species	Human
Gene ID	388531
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human R9AP/RGS9BP. (NP_997274.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q6ZS82
Protein Weight	25kDa
Shipping	Ice bag