


Anti-Rab10 Rabbit Polyclonal Antibody For Wb And Elisa - P61026

Item Number: CM0027732



Introduction

Rabbit polyclonal antibody against human RAB10, suitable for WB and ELISA. Cross-reacts with mouse. Target protein is a small GTPase involved in Golgi-to-plasma membrane transport, ciliogenesis, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	RAB10 Rabbit pAb
Remarks/Alias	RAB10
Species	Human
GeneID (Human)	10890
GeneID	10890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-198 of human RAB10 (NP_057215.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P61026
Protein Weight	23kDa
Shipping	Ice bag