


Anti-Rab11A Polyclonal Antibody For Wb, If/Icc, Elisa - P62491

Item Number: CM0015817



Introduction

Rabbit polyclonal antibody specific for human RAB11A, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with rat. Ideal for studying endocytic recycling and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	RAB11A Rabbit pAb
Remarks/Alias	YL8; RAB11A
Species	Human
Gene ID (Human)	8766
Gene ID	8766
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RAB11A (NP_004654.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P62491
Protein Weight	24kDa
Shipping	Ice bag