


Anti-Rab11A Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P62491

Item Number: CM0008703



Introduction

High-quality rabbit monoclonal antibody against human RAB11A (P62491), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying endocytic recycling and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	RAB11A Rabbit mAb
Remarks/Alias	YL8; RAB11A
Species	Human
Gene ID (Human)	8766
Gene ID	8766
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62491
Protein Weight	24kDa
Shipping	Ice bag