


Anti-Rab11A/B Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P62491 / Q15907

Item Number: CM0002335



Introduction

This rabbit monoclonal antibody against RAB11A/B is validated for WB, IF/ICC, IHC-P, and ELISA. It recognizes human, mouse, and rat RAB11A and RAB11B, key regulators of endocytic recycling and membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	RAB11A/B Rabbit mAb
Remarks/Alias	YL8; H-YPT3; NDAGSCW; RAB11A/B
Species	Human
GeneID (Human)	8766/9230
GeneID	8766/9230
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RAB11A/B (NP_004654.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62491/Q15907
Protein Weight	24KDa
Shipping	Ice bag