


Anti-Rab11A/Rab11B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62491 / Q15907

Item Number: CM0023083



Introduction

This rabbit polyclonal antibody detects RAB11A and RAB11B in human, mouse, and rat samples. Validated for WB, IHC-P, IF/ICC, and ELISA, it supports endocytic recycling and cell division research.

[Learn More](#)

Parameter	Value
Product Name	RAB11A/RAB11B Rabbit pAb
Species	Human
Gene ID (Human)	8766/9230
Gene ID	8766/9230
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RAB11A/RAB11B (NP_004654.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62491/Q15907
Shipping	Ice bag