


Anti-Rab11B Polyclonal Antibody For Wb, If/Icc, Elisa - Q15907

Item Number: CM0013127



Introduction

High-quality rabbit polyclonal antibody against human RAB11B (UniProt Q15907). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for endocytic recycling and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	RAB11B Rabbit pAb
Remarks/Alias	H-YPT3; NDAGSCW; RAB11B
Species	Human
GeneID (Human)	9230
GeneID	9230
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human RAB11B (NP_004209.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15907
Protein Weight	24kDa
Shipping	Ice bag