

Anti-Rab11Fip5 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bxf6

Item Number: CM0021045

Introduction

RAB11FIP5 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Recombinant immunogen (aa 354-653). Ideal for studying apical recycling endosome trafficking, insulin secretion, and V-ATPase transport.

[Learn More](#)

Parameter	Value
Product Name	RAB11FIP5 Rabbit pAb
Remarks/Alias	GAF1; pp75; RIP11; gaf-1; rab11-FIP5; RAB11FIP5
Species	Human
Gene ID	26056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 354-653 of human RAB11FIP5 (NP_056285.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXF6
Protein Weight	70kDa
Shipping	Ice bag