

Anti-Rab12 Monoclonal Antibody For Ip And Elisa - Q6Iq22

Item Number: CM0008184



Introduction

High-quality anti-Rab12 monoclonal antibody for immunoprecipitation (IP) and ELISA applications. Validated for human samples, targets the small GTPase Rab-12 involved in endosomal trafficking and autophagy. Ideal for IVD and research use.

[Learn More](#)

Parameter	Value
Product Name	Rab12 Rabbit mAb [MJFF-A-3F6]
Remarks/Alias	RAB12
Species	Human
Gene ID (Human)	201475
Gene ID	201475
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-244 of human Rab12 protein (Q6IQ22).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IP, ELISA
Cross Reactivity	Human
SWISS	Q6IQ22
Protein Weight	27kDa
Shipping	Ice bag