

Anti-Rab13 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P51153

Item Number: CM0014718

Introduction

Rabbit polyclonal antibody against human RAB13, validated for WB, IF/ICC, ELISA. Ideal for studying tight junction assembly, endocytic recycling, and GLUT4 translocation.

[Learn More](#)

Parameter	Value
Product Name	RAB13 Rabbit pAb
Remarks/Alias	GIG4; RAB13
Species	Human
GeneID (human)	5872
GeneID	5872
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P51153
Protein Weight	23kDa
Shipping	Ice bag