

Anti-Rab14 Rabbit Polyclonal Antibody For Wb And Elisa - P61106

Item Number: CM0014224



Introduction

Anti-RAB14 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human and mouse RAB14, a Rab small GTPase regulating Golgi-endosome traffic, endosomal tethering, and cell adhesion. Suitable for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	RAB14 Rabbit pAb
Remarks/Alias	FBP, RAB-14, RAB14
Species	Human
GeneID (Human)	51552
GeneID	51552
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 39-140 of human RAB14 (NP_057406.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61106
Protein Weight	24kDa
Shipping	Ice bag