

Anti-Rab18 Polyclonal Antibody For Wb And Elisa - Q9Np72

Item Number: CM0030445



Introduction

Rabbit polyclonal antibody against human RAB18, validated for WB and ELISA. Detects endogenous RAB18 in human, mouse, and rat samples. Ideal for studying intracellular membrane trafficking, ER-lipid droplet contacts, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	RAB18 Rabbit pAb
Remarks/Alias	WARBM3; RAB18LI1; RAB18
Species	Human
Gene ID (Human)	22931
Gene ID	22931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RAB18 (NP_067075.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NP72
Protein Weight	23kDa
Shipping	Ice bag