

Anti-Rab1A Polyclonal Antibody For Wb, If/Icc, Elisa - P62820

Item Number: CM0019195



Introduction

Rabbit polyclonal antibody against human RAB1A, suitable for WB, IF/ICC, ELISA; cross-reactive with mouse and rat. Ideal for studying ER-Golgi trafficking, autophagy, and cell migration. 23 kDa.

[Learn More](#)

Parameter	Value
Product Name	RAB1A Rabbit pAb
Remarks/Alias	RAB1; YPT1; RAB1A
Species	Human
Gene ID (Human)	5861
Gene ID	5861
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-205 of human RAB1A (NP_004152.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62820
Protein Weight	23kDa
Shipping	Ice bag