

Anti-Rab1A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62820

Item Number: CM0013224



Introduction

Rabbit polyclonal antibody targeting human RAB1A (NP_004152.1). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for ER-to-Golgi trafficking and autophagy studies.

[Learn More](#)

Parameter	Value
Product Name	RAB1A Rabbit pAb
Remarks/Alias	RAB1; YPT1; RAB1A
Species	Human
Gene ID (Human)	5861
Gene ID	5861
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human RAB1A (NP_004152.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P62820
Protein Weight	23kDa
Shipping	Ice bag