


Anti-Rab1B Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0U4

Item Number: CM0027961



Introduction

High-quality rabbit polyclonal antibody against human RAB1B (NP_112243.1), validated for Western blot and ELISA. Suitable for detecting RAB1B in human and mouse samples, with a molecular weight of 22 kDa. Ideal for studying intracellular membrane trafficking and autophagy.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RAB1B Rabbit pAb
Remarks/Alias	RAB1B
Species	Human
Gene ID (Human)	81876
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 7-201 of human RAB1B (NP_112243.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H0U4
Protein Weight	22kDa
Shipping	Ice bag