


Anti-Rab1B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H0U4

Item Number: CM0022199



Introduction

Rabbit polyclonal antibody against human RAB1B (ras-related protein Rab-1B). Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Suitable for studying ER-Golgi trafficking, autophagy, and GTPase signaling.

[Learn More](#)

Specification	Details
Product Name	RAB1B Rabbit pAb
Remarks/Alias	RAB1B; ras-related protein Rab-1B
Species	Human
GeneID (Human)	81876
GeneID	81876
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 7-201 of human RAB1B (NP_112243.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0U4
Protein Weight	22kDa
Shipping	Ice bag