

Anti-Rab20 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nx57

Item Number: CM0019379

Introduction

Rabbit polyclonal antibody targeting human RAB20. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Supports research into membrane trafficking, phagosome maturation, and pancreatic cancer biology.

[Learn More](#)

Parameter	Value
Product Name	RAB20 Rabbit pAb
Remarks/Alias	RAB20
Species	Human
Gene ID (Human)	55647
Gene ID	55647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-190 of human RAB20 (NP_060287.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q9NX57
Protein Weight	26kDa
Shipping	Ice bag