

Anti-Rab21 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9UI25

Item Number: CM0015262



Introduction

Rabbit polyclonal antibody against human RAB21 (Ras-related protein Rab-21), validated for WB, IF/ICC, and ELISA. Detects endogenous RAB21 in human, mouse, and rat samples at ~24 kDa. Ideal for studies of membrane trafficking, integrin recycling, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	RAB21 Rabbit pAb
Remarks/Alias	RAB21
Species	Human
GeneID (Human)	23011
GeneID	23011
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL25
Protein Weight	24kDa
Shipping	Ice bag