

Anti-Rab22A Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9UI26

Item Number: CM0006818

Introduction

Rabbit monoclonal anti-RAB22A antibody validated for WB, IP, ELISA in Human, Mouse, Rat. Recombinant immunogen (aa 1-100). Targets ras-related protein Rab-22A, a key endocytic trafficking GTPase.

[Learn More](#)

Parameter	Value
Product Name	RAB22A Rabbit mAb
Remarks/Alias	RAB22A; ras-related protein Rab-22A
Species	Human
GeneID (Human)	57403
GeneID	57403
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human RAB22A (NP_065724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL26
Protein Weight	22 kDa
Shipping	Ice bag