


Anti-Rab23 Rabbit Polyclonal Antibody For Wb, Elisa - Q9ULc3

Item Number: CM0027637



Introduction

Rabbit polyclonal antibody against human RAB23 for Western blot and ELISA. Validated in human, mouse, and rat. Targets a 27kDa protein involved in ciliary trafficking and hedgehog signaling.

[Learn More](#)

Parameter	Value
Product Name	RAB23 Rabbit pAb
Remarks/Alias	HSPC137; RAB23
Species	Human
GeneID (Human)	51715
GeneID	51715
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULC3
Protein Weight	27kDa
Shipping	Ice bag