

Anti-Rab24 Polyclonal Antibody For Wb, Ip, Elisa - Q969Q5

Item Number: CM0026655

Introduction

Anti-RAB24 Rabbit Polyclonal Antibody validated for WB, IP, ELISA in human, mouse, rat. Targets Ras-related protein Rab-24 (UniProt Q969Q5), involved in autophagy and intracellular membrane trafficking. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	RAB24 Rabbit pAb
Remarks/Alias	RAB24
Species	Human
GeneID (Human)	53917
GeneID	53917
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-203 of human RAB24 (NP_001026847.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969Q5
Protein Weight	23kDa
Shipping	Ice bag