

## Anti-Rab27B Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O00194

Item Number: CM0012471



### Introduction

Rabbit monoclonal antibody against human RAB27B, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying late endocytic pathway regulation and intracellular trafficking.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | RAB27B Rabbit mAb                                                                                                 |
| Remarks/Alias    | C25KG                                                                                                             |
| Species          | Human                                                                                                             |
| GeneID (human)   | 5874                                                                                                              |
| GeneID           | 5874                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human RAB27B (NP_004154.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O00194                                                                                                            |
| Protein Weight   | 25kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |