


Anti-Rab28 Rabbit Polyclonal Antibody For Wb, Elisa - P51157

Item Number: CM0022222



Introduction

RAB28 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human RAB28, a small GTPase involved in membrane trafficking and cone cell outer segment shedding. Suitable for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	RAB28 Rabbit pAb
Remarks/Alias	CORD18; RAB28
Species	Human
GeneID (Human)	9364
GeneID	9364
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51157
Protein Weight	25kDa
Shipping	Ice bag