


## Anti-Rab2A Polyclonal Antibody For Wb, Elisa - P61019

Item Number: CM0020419



### Introduction

Rabbit polyclonal antibody targeting human RAB2A (Ras-related protein Rab-2A). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Regulates ER-Golgi transport and autophagy. Research use only.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | RAB2A Rabbit pAb                                                                                                                      |
| Remarks/Alias    | LHX; RAB2; RAB2A                                                                                                                      |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 5862                                                                                                                                  |
| Gene ID          | 5862                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human RAB2A (NP_002856.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | P61019                                                                                                                                |
| Protein Weight   | 24kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |