


Anti-Rab2A Rabbit Polyclonal Antibody For Wb And Elisa - P61019

Item Number: CM0016208



Introduction

High-quality rabbit polyclonal antibody targeting human RAB2A for Western blot and ELISA. Validated in human, mouse, and rat samples. Ideal for autophagy and intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	RAB2A Rabbit pAb
Remarks/Alias	LHX, RAB2, RAB2A
Species	Human
GeneID (Human)	5862
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human RAB2A (NP_002856.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61019
Protein Weight	24kDa
Shipping	Ice bag