

Anti-Rab31 Polyclonal Antibody For Wb, Elisa - Q13636

Item Number: CM0030431



Introduction

Rabbit polyclonal antibody targeting human RAB31 (Rab22B) for western blot and ELISA. Cross-reacts with mouse and rat. Validated with recombinant immunogen. Explore the role of this Rab GTPase in Golgi function and endocytosis.

[Learn More](#)

Parameter	Value
Product Name	RAB31 Rabbit pAb
Remarks/Alias	Rab22B; RAB31
Species	Human
Gene ID (Human)	11031
Gene ID	11031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human RAB31 (NP_006859.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13636
Protein Weight	22kDa
Shipping	Ice bag