

Anti-Rab31 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13636

Item Number: CM0026549



Introduction

RAB31 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Explore the role of RAB31 GTPase in Golgi function, endocytosis, and phagosome maturation.

[Learn More](#)

Parameter	Value
Product Name	RAB31 Rabbit pAb
Remarks/Alias	Rab22B; RAB31
Species	Human
GeneID (Human)	11031
GeneID	11031
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human RAB31 (NP_006859.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13636
Protein Weight	22kDa
Shipping	Ice bag