

Anti-Rab33A (Rabs10) Rabbit Polyclonal Antibody For Wb And Elisa - Q14088

Item Number: CM0020771

Introduction

Rabbit polyclonal antibody against human RAB33A, validated for Western blot and ELISA with cross-reactivity to mouse and rat. Ideal for studying autophagosome formation and membrane trafficking in lymphoid cells.

[Learn More](#)

Parameter	Value
Product Name	RAB33A Rabbit pAb
Remarks/Alias	RabS10; RAB33A
Species	Human
GeneID (Human)	9363
GeneID	9363
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-237 of human RAB33A (NP_004785.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14088
Protein Weight	27kDa
Shipping	Ice bag