

Anti-Rab3Gap1 Polyclonal Antibody For Wb, Elisa - Q15042

Item Number: CM0018854



Introduction

RAB3GAP1 Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, rat. The catalytic subunit of the Rab3GAP complex with roles in neurodevelopment and intracellular trafficking.

[Learn More](#)

Parameter	Value
Product Name	RAB3GAP1 Rabbit pAb
Remarks/Alias	P130; MARTS2; WARBM1; RAB3GAP; RAB3GAP130; RAB3GAP1
Species	Human
GeneID (Human)	22930
GeneID	22930
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-780 of human RAB3GAP1 (NP_001165906.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15042
Protein Weight	111kDa
Shipping	Ice bag