


Anti-Rab3Ip Polyclonal Antibody For Wb, Elisa - Q96Qf0

Item Number: CM0030381



Introduction

RAB3IP Rabbit pAb validated for Western blot and ELISA in human, mouse, and rat. Target RAB3IP is a GEF for RAB8A/B, involved in vesicle trafficking and ciliogenesis. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	RAB3IP Rabbit pAb
Remarks/Alias	RABIN3; RABIN8; RAB3IP
Species	Human
GeneID (Human)	117177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human RAB3IP (NP_071901.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96QF0
Protein Weight	53kDa
Shipping	Ice bag