


Anti-Rab45 Rabbit Monoclonal Antibody For Wb And Elisa - Q8Iz41

Item Number: CM0006777



Introduction

CamelBio's RAB45 rabbit monoclonal antibody (clone RAB45, immunogen: N-terminal 1-120 aa) is validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for studying dynein-mediated transport, endosomal trafficking, and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	RAB45 Rabbit mAb
Remarks/Alias	TSG; RAB45
Species	Human
Gene ID (Human)	158158
Gene ID	158158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human RAB45 (NP_689786.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IZ41
Protein Weight	83kDa
Shipping	Ice bag