


Anti-Rab4A Monoclonal Antibody For Wb - P20338

Item Number: CM0010825



Introduction

Rabbit monoclonal antibody targeting human RAB4A (P20338). Validated for WB and ELISA, with cross-reactivity to mouse. Ideal for studying intracellular trafficking, VEGFR2 signaling, and endosomal pathways.

[Learn More](#)

Parameter	Details
Product Name	RAB4A Rabbit mAb
Remarks/Alias	RAB4; HRES1; HRES-1; HRES-1/RAB4; RAB4A
Species	Human
GeneID (Human)	5867
GeneID	5867
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20338
Protein Weight	24kDa
Shipping	Ice bag