

Anti-Rab5A Monoclonal Antibody For Wb, Ihc-P, Elisa - P20339

Item Number: CM0007442



Introduction

High-quality rabbit monoclonal antibody targeting human RAB5A, a key regulator of endosomal trafficking. Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for endocytic trafficking and membrane fusion research.

[Learn More](#)

Parameter	Value
Product Name	RAB5A Rabbit mAb
Remarks/Alias	RAB5A; RAB5; ras-related protein Rab-5A
Species	Human
Gene ID (Human)	5868
Gene ID	5868
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 115-215 of human RAB5A (NP_004153.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20339
Protein Weight	22kDa/23kDa
Shipping	Ice bag