

Anti-Rab5A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P20339

Item Number: CM0013349



Introduction

RAB5A Rabbit Polyclonal Antibody (P20339) validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets early endosomal trafficking regulator RAB5A; ideal for membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	RAB5A Rabbit pAb
Remarks/Alias	RAB5; RAB5A
Species	Human
GeneID (Human)	5868
GeneID	5868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human RAB5A (NP_004153.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20339
Protein Weight	24kDa
Shipping	Ice bag