

Anti-Rab5B Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P61020

Item Number: CM0026021



Introduction

RAB5B Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, IF-P, and ELISA in human, mouse, and rat. Targets Ras-related protein Rab-5B, a key regulator of early endocytic trafficking for membrane dynamics research.

[Learn More](#)

Parameter	Value
Product Name	RAB5B Rabbit pAb
Remarks/Alias	RAB5B
Species	Human
GeneID (Human)	5869
GeneID	5869
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RAB5B (NP_002859.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61020
Protein Weight	24kDa
Shipping	Ice bag