

Anti-Rab5C Polyclonal Antibody For Wb, Elisa - P51148

Item Number: CM0022904



Introduction

Rabbit polyclonal antibody against human RAB5C (SWISS P51148) for WB and ELISA applications. Immunogen: recombinant fusion protein corresponding to amino acids 1-216. Cross-reacts with human; protein weight 23kDa. Ideal for endocytic trafficking research.

[Learn More](#)

Parameter	Value
Product Name	RAB5C Rabbit pAb
Remarks/Alias	RABL; L1880; RAB5L; RAB5CL; RAB5C
Species	Human
Gene ID (Human)	5878
Gene ID	5878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-216 of human RAB5C (NP_958842.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51148
Protein Weight	23kDa
Shipping	Ice bag