

Anti-Rab5C / Rabl Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P51148

Item Number: CM0003225



Introduction

Rabbit monoclonal antibody targeting human RAB5C/RABL, suitable for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Molecular weight: 23kDa. Ideal for studying early endocytic trafficking.

[Learn More](#)

Parameter	Value
Product Name	RAB5C / RABL Rabbit mAb
Remarks/Alias	RABL; L1880; RAB5L; RAB5CL; RAB5C / RABL
Species	Human
GeneID (Human)	5878
GeneID	5878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-216 of human RAB5C / RABL (P51148)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51148
Protein Weight	23kDa
Shipping	Ice bag