

Anti-Rab7 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P51149

Item Number: CM0005828



Introduction

Mouse monoclonal antibody targeting human RAB7, suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Essential for studying endo-lysosomal trafficking and related diseases.

[Learn More](#)

Parameter	Value
Product Name	RAB7 Mouse mAb
Remarks/Alias	RAB7; CMT2B; PRO2706
Species	Human
Gene ID (Human)	7879
Gene ID	7879
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 108-207 of human RAB7 (NP_004628.4).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51149
Protein Weight	23 kDa
Shipping	Ice bag