

Anti-Rab7 Polyclonal Antibody For Wb, Ip, Elisa - P51149

Item Number: CM0013815



Introduction

High-quality rabbit polyclonal antibody against RAB7 (Ras-related protein Rab-7a) validated for WB, IP, ELISA in human, mouse, rat. Ideal for endosomal trafficking, autophagy, and vesicular transport research.

[Learn More](#)

Parameter	Value
Product Name	RAB7 Rabbit pAb
Remarks/Alias	RAB7; CMT2B; PRO2706
Species	Human
GeneID (Human)	7879
GeneID	7879
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-207 of human RAB7 (NP_004628.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51149
Protein Weight	23kDa
Shipping	Ice bag