

Anti-Rab7 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P51149

Item Number: CM0003289



Introduction

Recombinant rabbit monoclonal antibody targeting human RAB7. Validated for WB, IP, IF/ICC, IHC-P, ELISA; reacts with human, mouse, rat. Essential for endosomal trafficking studies; suitable for IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	RAB7 Rabbit mAb
Remarks/Alias	RAB7; CMT2B; PRO2706
Species	Human
GeneID (Human)	7879
GeneID	7879
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 108-207 of human RAB7 (NP_004628.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51149
Protein Weight	23kDa
Shipping	Ice bag