


Anti-Rab8A Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P61006

Item Number: CM0002670



Introduction

Rabbit monoclonal antibody targeting human RAB8A (MEL, RAB8). Validated for WB, IF/ICC, IP, ELISA; cross-reacts with mouse and rat. Suitable for intracellular trafficking, ciliogenesis, and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	RAB8A Rabbit mAb
Remarks/Alias	MEL; RAB8; RAB8A
Species	Human
GeneID (Human)	4218
GeneID	4218
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-207 of human RAB8A (P61006).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61006
Protein Weight	24kDa
Shipping	Ice bag