


Anti-Rab8A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61006

Item Number: CM0019214



Introduction

RAB8A Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat RAB8A (24 kDa). Ideal for studying intracellular trafficking, ciliogenesis, autophagy, and epithelial polarization.

[Learn More](#)

Parameter	Value
Product Name	RAB8A Rabbit pAb
Remarks/Alias	MEL; RAB8; RAB8A
Species	Human
GeneID (Human)	4218
GeneID	4218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-207 of human RAB8A (NP_005361.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61006
Protein Weight	24kDa
Shipping	Ice bag