

Anti-Rab9A Monoclonal Antibody For Wb, If/Icc, Elisa - P51151

Item Number: CM0009118



Introduction

Rabbit monoclonal antibody targeting human RAB9A, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. RAB9A is a key regulator of endosome-to-TGN trafficking. Ideal for intracellular transport studies.

[Learn More](#)

Parameter	Value
Product Name	RAB9A Rabbit mAb
Remarks/Alias	RAB9; RAB9A
Species	Human
GeneID (Human)	9367
GeneID	9367
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human RAB9A (NP_004242.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51151
Protein Weight	23kDa
Shipping	Ice bag