

Anti-Rab9P40 Monoclonal Antibody For Wb, Elisa - Q7Z6M1

Item Number: CM0003464

Introduction

Recombinant rabbit monoclonal antibody specific for human RAB9P40 (Rab9 effector p40), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying endosome-to-TGN transport in diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	RAB9P40 Rabbit mAb
Remarks/Alias	p40; RAB9P40; bA65N13.1
Species	Human
Gene ID	10244
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-372 of human RAB9P40 (Q7Z6M1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z6M1
Protein Weight	41kDa
Shipping	Ice bag