

## Anti-Rabac1 Polyclonal Antibody For Wb, Elisa - Q9Ui14

Item Number: CM0019494

### Introduction

Rabbit polyclonal antibody targeting human RABAC1 (PRA1), suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Useful for studies of Golgi vesicle trafficking and Rab protein regulation.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | RABAC1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | PRA1; YIP3; PRAF1; RABAC1                                                                                          |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 10567                                                                                                              |
| Gene ID          | 10567                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human RABAC1 (NP_006414.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q9Ui14                                                                                                             |
| Protein Weight   | 21kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |