


Anti-Rabgap1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y3P9

Item Number: CM0017989



Introduction

Rabbit polyclonal antibody targeting RABGAP1, a Rab GTPase-activating protein involved in centrosome function and cell cycle progression. Suitable for WB, IHC-P, and ELISA; reacts with human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	RABGAP1 Rabbit pAb
Remarks/Alias	GAPCENA; TBC1D11; RABGAP1
Species	Human
GeneID (Human)	23637
GeneID	23637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 830-1069 of human RABGAP1 (NP_036329.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y3P9
Protein Weight	122kDa
Shipping	Ice bag